# Demotec 95

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**Demotec 95** is a specialised resin specifically developed for the adjuvant treatment of lameness in cattle and represents a significant advance on Demotec 90 and other old fashioned bonding agents available.

**Demotec 95** possesses the ideal characteristics for this purpose as the putty-like mass formed a few seconds after mixing, is malleable yet not sticky, so does not adhere to the hands of the operator.

**Demotec 95** will bond a wooden bloc onto the sole of the sound claw, so that it transfers the weight onto the good claw, allowing the affected claw to rest. Additionally it will be raised out of the dirt, thus making medical and surgical treatments more effective. This consequently relieves the animal of pain and directly reduces the degree of lameness.

**Demotec 95** also is available with large wooden blocks.

#### **Preparation of the hoof**

Should it be necessary, the good claw should be trimmed and cleaned of dirt and loose particles using a hoof knife or angle cutter. It is recommended that the sole of the claw should then be cleaned with spirit or alcohol in order to improve the bonding of the resin onto the claw. Having prepared the hoof, the mixed resin should then be applied to both the claw and the wooden bloc.

#### **Preparation of Demotec 95**

This specialised resin consists of a powder and liquid which should be mixed in a ratio of 2 : 1. For ease of mixing the powder is presented in sachets of 70 g, whilst a beaker graduated to 35 mls is provided for the liquid.

- Measure 35 mls of liquid and pour into the beaker
- Empty the contents of one sachet of powder (70 g) into the liquid
- Mix thoroughly with the wooden spatula provided.

#### **Special note**

As the material is not liquid as e.g. Demotec 90 or other old fashioned products, the whole operation should be undertaken without any mess. The mixing beaker is intended for use on several occasions; normally the mixed resin will come away from the beaker, but when set any remaining will separate from the beaker with slight pressure.

#### **Processing and application of Demotec 95**

After a few seconds an homogenous putty-like mass will be formed which will conveniently adhere to the spatula. This kneadable mass can then be modelled very easily for about 1 – 2 minutes. After 4 – 5 minutes the resin will have set completely and when cold it can withstand very heavy loads.

It is recommended that the sole of the claw be covered with a thin layer of resin as this will optimise bonding; then the resin should be applied, at the thickness required to the surface of the wooden bloc and once in position, modelled into any desired shape, wrapping the resin round the sides of the claw and the block and into the circumference grooves of the block.

#### **Removing the wooden bloc**

When the affected claw is better, the wooden bloc can be removed using a chisel or another cutting tool.

#### **Setting-times – affected by temperature**

Because on occasions it is necessary to treat animals outdoors it should be noted that the setting time of the resin is significantly affected by temperature.

Times quoted in this document refer to an ambient temperature of 21°C. The resin may be heated in the summer or cooled in the winter. Lower temperatures result in a longer setting time whilst higher temperatures a shorter setting time. To overcome this problem, we have developed the Demotec Accelerator for using during the winter which significantly reduces the polymerisation time.

temp.	modelling times	setting times	setting times with Demotec Accelerator
21 °C	1-2 min.	≈ 4 min.	≈ 2 min.
15 °C	2-3 min.	≈ 8 min.	≈ 4 min.
5 °C	4-6 min.	≈ 10 min.	≈ 6 min.

#### **Other applications**

Demotec 95, besides being used in the treatment of lameness in cattle has also been used for splinting fractures in small animals – because it can so easily be moulded and controlled.

#### **General**

Demotec 95 should be:

**Stored in a cool place**

**Protected from light**

**Kept out of the reach of children**

The warnings on the bottles should be noted and observed. Further details and DIN safety specifications are available on request.

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